

# INJECTOR SUCTION



PART No. 4425



PART No. 4426

The Therapy Equipment Injector Suction Assemblies are ideal for all types of emergency suction use, and can be fitted to all types of trolleys.

The Injector Suction Unit is designed to provide maximum vacuum with the minimum usage of gas.

The Unit will provide: Vacuum - up to 400mm HG/55 Kpa.

Gas Usage - 19 LPM.

In line with the latest British and ISO Standards (BS7259/ISO 10079) the Unit has a Pressure Relief Valve to prevent any direct pressure being applied to the patient.

The Suction Connection is through the Pipeline Protector at the bottom of the Unit, the same as for Pipeline Controllers.

The Injector Suction is therefore protected from contamination.

The Suction Control Valve is situated at the rear of the Suction Controller Bonnet.

The Injector Suction is available in the following configurations:

Part Number 4420 - Injector Suction and Flowmeter on Backbar.

Part Number 4425 - Injector Suction/Flowmeter Rail Mounted.

Part Number 4426 - Injector Suction/Flowgauge Rail Mounted.

Part Number 4427 - As for 4425 with Reg and Schrader.

Part Number 4428 - As for 4426 with Reg and Schrader.


Part Number 4508 - 1/2 Litre Jar Assembly.

## **THErapy EQUIPMENT LTD**

UNIT 1, CRANBOURNE AVENUE  
CRANBOURNE INDUSTRIAL ESTATE  
POTTERS BAR, HERTS. EN6 3JN

TEL: 01707 652270 FAX: 01707 652622

E-mail: [sales@therapiequipment.co.uk](mailto:sales@therapiequipment.co.uk)  
Website: [www.therapiequipment.co.uk](http://www.therapiequipment.co.uk)



0086



FM 37970